

Notch1 Rabbit mAb


ABZ4752

Product Overview

Name	Notch1 Rabbit mAb
Accession(Primary)	P46531
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 273 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:100-1:200

Target information (P46531)

Gene Name	NOTCH1
About Gene Symbol	NOTCH1 (Notch Receptor 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Notch1 in C6, HeLa lysates using Notch1 antibody.