

LRP6 Rabbit mAb


ABZ4754

Product Overview

Name	LRP6 Rabbit mAb
Accession(Primary)	O75581
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 180 kDa; Observed MW: 180-230 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O75581)

Gene Name	LRP6
About Gene Symbol	



Western blot analysis of LRP6 in HeLa and HepG2
lysates using LRP6 antibody.

Copyright© 2024 Zen-bioscience Limited.