


[KO]beta Dystroglycan Rabbit mAb

ABZ4757

Product Overview

Name	[KO]beta Dystroglycan Rabbit mAb
Accession(Primary)	Q14118
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 97 kDa; Observed MW: 42-44 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q14118)

Gene Name DAG1**About Gene Symbol**

Western blot analysis of beta Dystroglycan in MCF-7, HepG2, Lncap lysates using beta Dystroglycan antibody.

