

[KD]KDM2A Rabbit mAb


ABZ4758

Product Overview

Name	[KD]KDM2A Rabbit mAb
Accession(Primary)	Q9Y2K7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 133 kDa; Observed MW: 133 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:500-1:1000; IP-1:50-1:100

Target information (Q9Y2K7)

Gene Name	KDM2A
About Gene Symbol	KDM2A (Lysine Demethylase 2A) is a transcription factor involved in chromatin remodeling and epigenetics. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and id



Western blot analysis of FBXL11 in Jurkat lysates
using KDM2A antibody.