

[KO]IL-1RAcP Rabbit mAb

ABZ4760


Product Overview

Name	[KO]IL-1RAcP Rabbit mAb
Accession(Primary)	Q9NPH3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 65 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:50-1:100; IP-1:20-1:50

Target information (Q9NPH3)

Gene Name IL1RAP

About Gene Symbol



Western blot analysis of IL-1RAcP in HL-60, U251, U87-MG lysates using IL-1 Receptor Accessory Protein antibody.

