

[KD]PAIP1 Rabbit mAb


ABZ4761

Product Overview

Name	[KD]PAIP1 Rabbit mAb
Accession(Primary)	Q9H074
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; FC-1:10-1:100; IP-1:10-1:100

Target information (Q9H074)

Gene Name	PAIP1
About Gene Symbol	PAIP1 (Poly(A) Binding Protein Interacting Protein 1) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PAIP1 in HeLa, 293T, mouse brain lysates using PAIP1 antibody.