

[KO]LPCAT1 Rabbit mAb


ABZ4762

Product Overview

Name	[KO]LPCAT1 Rabbit mAb
Accession(Primary)	Q8NF37
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20-1:50

Target information (Q8NF37)

Gene Name	LPCAT1
About Gene Symbol	LPCAT1 (Lysophosphatidylcholine Acyltransferase 1) is a mitochondrial acyltransferase enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate



Western blot analysis of LPCAT1 in A549 lysates
using LPCAT1 antibody.