

[KD]TORC2 Rabbit mAb


ABZ4763

Product Overview


Name	[KD]TORC2 Rabbit mAb
Accession(Primary)	Q53ET0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (Q53ET0)

Gene Name	CRTC2
About Gene Symbol	CRTC2 (CREB Regulated Transcription Coactivator 2) is a transcriptional coactivator involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate c



Western blot analysis of TORC2 in A549, HL-60, U251 lysates using CREB Regulated Transcription Coactivator 2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using TORC2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.