

[KO]CPSF6 Rabbit mAb


ABZ4764

Product Overview

Name	[KO]CPSF6 Rabbit mAb
Accession(Primary)	Q16630
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 59 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:10-1:100

Target information (Q16630)

Gene Name	CPSF6
About Gene Symbol	CPSF6 (Cleavage And Polyadenylation Specific Factor 6) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to



Western blot analysis of CPSF6 in rat Brain, MCF-7, 3T3, C6 lysates using CPSF6 antibody.