

[KO]EFEMP1 Rabbit mAb

ABZ4765

Product Overview

Name	[KO]EFEMP1 Rabbit mAb
Accession(Primary)	Q12805
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50

Target information (Q12805)

Gene Name	EFEMP1
About Gene Symbol	EFEMP1 (EGF Containing Fibulin Extracellular Matrix Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

