

[KO]DAB2 Rabbit mAb


ABZ4766

Product Overview


Name	[KO]DAB2 Rabbit mAb
Accession(Primary)	P98082
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 82 kDa; Observed MW: 96 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IP-1:20-1:50

Target information (P98082)

Gene Name	DAB2
About Gene Symbol	DAB2 (DAB Adaptor Protein 2) is a signaling adaptor protein involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of DAB2 in HepG2, U2OS, U251 lysates using DAB2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using DAB2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.