

[KD]TCTEX-1 Rabbit mAb


ABZ4767

Product Overview


Name	[KD]TCTEX-1 Rabbit mAb
Accession(Primary)	P63172
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (P63172)


Gene Name	DYNLT1
About Gene Symbol	DYNLT1 (Dynein Light Chain Tctex-Type 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Immunocytochemistry analysis of Dynein Light Chain Tctex Type 1 (green) in Hela using Dynein Light Chain Tctex Type 1 antibody, and DAPI (blue).



Western blot analysis of TCTEX1 in 3T3, C6 lysates using TCTEX1 antibody.



Western blot analysis of TCTEX1 in HeLa, A549, HL-60 lysates using TCTEX1 antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TCTEX1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.