

[KD]UFM1 Rabbit mAb


ABZ4768

Product Overview


Name	[KD]UFM1 Rabbit mAb
Accession(Primary)	P61960
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 9 kDa; Observed MW: 9 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; IP-1:10-1:50

Target information (P61960)


Gene Name	UFM1
About Gene Symbol	UFM1 (Ubiquitin Fold Modifier 1) is a metabolic enzyme involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of UFM1 in 3T3, Raw264.7 lysates using UFM1 antibody.



Western blot analysis of UFM1 in Jurkat lysates using UFM1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using UFM1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.