


# [KO]HspA2 Rabbit mAb

ABZ4772


## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | [KO]HspA2 Rabbit mAb                       |
| <b>Accession(Primary)</b>   | P54652                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 70 kDa; Observed MW: 70 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:2000-1:5000      |


## Target information (P54652)

**Gene Name** HSPA2**About Gene Symbol**

Western blot analysis of HSPA2 in A549, Hacat lysates using HspA2 antibody.



Western blot analysis of HSPA2 in rat Brain, C6, 3T3 lysates using HSPA2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using HSPA2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.