

[KD]D4 GDI Rabbit mAb


ABZ4773

Product Overview

Name	[KD]D4 GDI Rabbit mAb
Accession(Primary)	P52566
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (P52566)

Gene Name	ARHGDIB
About Gene Symbol	



Western blot analysis of D4 GDI in Jurkat, U937 lysates using D4 GDI antibody.