

[KD]CENPA Rabbit mAb


ABZ4774

Product Overview

Name	[KD]CENPA Rabbit mAb
Accession(Primary)	P49450
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P49450)

Gene Name	CENPA
About Gene Symbol	CENPA (Centromere Protein A) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics and adaptive immunity. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CENPA in K562, HeLa
lysates using CENPA antibody.