


[KD]NMDAR2A Rabbit mAb

ABZ4779

Product Overview

Name	[KD]NMDAR2A Rabbit mAb
Accession(Primary)	P35436
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 165 kDa; Observed MW: 165 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P35436)

Gene Name Grin2A**About Gene Symbol**

Western blot analysis of NMDAR2A in mouse lung lysates using NMDAR2A antibody.

