

[KO]ARPC5/p16 ARC Rabbit mAb


ABZ4783

Product Overview


Name	[KO]ARPC5/p16 ARC Rabbit mAb
Accession(Primary)	O15511
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; IP-1:20-1:100

Target information (O15511)

Gene Name	ARPC5
About Gene Symbol	



Western blot analysis of p16 ARC in rat Brain, C6, 3T3, Hela lysates using p16 ARC antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using p16 ARC antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.