

[KD]CD168 Rabbit mAb


ABZ4787

Product Overview

Name	[KD]CD168 Rabbit mAb
Accession(Primary)	O75330
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:10-1:100

Target information (O75330)

Gene Name	HMMR
About Gene Symbol	HMMR (Hyaluronan Mediated Motility Receptor) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CD168 in K562 lysates
using CD168 antibody.