

[KD]Cyclin T1 Rabbit mAb


ABZ4790

Product Overview


Name	[KD]Cyclin T1 Rabbit mAb
Accession(Primary)	O60563
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (O60563)

Gene Name	CCNT1
About Gene Symbol	CCNT1 (Cyclin T1) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cyclin T1 in mouse brain, rat brain lysates using Cyclin T1 antibody.



Western blot analysis of Cyclin T1 in HeLa, Jurkat, K562 lysates using Cyclin T1 antibody.