

[KO]E1B AP5 Rabbit mAb


ABZ4791

Product Overview

Name	[KO]E1B AP5 Rabbit mAb
Accession(Primary)	Q9BUJ2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 96 kDa; Observed MW: 96 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500

Target information (Q9BUJ2)

Gene Name	HNRNPUL1
About Gene Symbol	HNRNPUL1 (Heterogeneous Nuclear Ribonucleoprotein U Like 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of E1BAP5 in K562, C6
lysates using E1B AP5 antibody.