

[KO]EMAP II Rabbit mAb


ABZ4792

Product Overview


Name	[KO]EMAP II Rabbit mAb
Accession(Primary)	Q12904
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q12904)

Gene Name	AIMP1
About Gene Symbol	AIMP1 (Aminoacyl TRNA Synthetase Complex Interacting Multifunctional Protein 1) is a biosynthetic enzyme involved in JAK-STAT/cytokine signaling and angiogenesis. It is an important research target in immunology, inflammation, and cancer research and is w



Western blot analysis of EMAP II/AIMP1 in HT-1080 lysates using EMAP II antibody.



Western blot analysis of EMAP II in C6 lysates using EMAP II antibody.