

## [KD]HNRNPA1L2 Rabbit mAb


ABZ4794

### Product Overview

|                             |                                                                     |
|-----------------------------|---------------------------------------------------------------------|
| <b>Name</b>                 | [KD]HNRNPA1L2 Rabbit mAb                                            |
| <b>Accession(Primary)</b>   | Q32P51                                                              |
| <b>Host Species</b>         | Rabbit                                                              |
| <b>Clonality</b>            | Monoclonal                                                          |
| <b>Isotype</b>              | IgG                                                                 |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                   |
| <b>Calculated MW</b>        | Calculated MW: 34 kDa; Observed MW: 34 kDa                          |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:100-1:200; FC-1:20 |

### Target information (Q32P51)

|                          |                                                                                                                                                                                                                                                                                                                                                |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | HNRNPA1L2                                                                                                                                                                                                                                                                                                                                      |
| <b>About Gene Symbol</b> | HNRNPA1L2 (Heterogeneous Nuclear Ribonucleoprotein A1 Like 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Heterogeneous Nuclear Ribonucleoprotein (A1like) in K562, C6, 3T3, HeLa lysates using HNRNPA1L2 antibody.