

[KD]PDLIM7 Rabbit mAb


ABZ4797

Product Overview

Name	[KD]PDLIM7 Rabbit mAb
Accession(Primary)	Q9NR12
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; FC-1:50-1:100

Target information (Q9NR12)

Gene Name	PDLIM7
About Gene Symbol	PDLIM7 (PDZ And LIM Domain 7) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PDLIM7 in HeLa lysates
using PDLIM7 antibody.