

[KO]E1B AP5 Rabbit mAb


ABZ4801

Product Overview



Name	[KO]E1B AP5 Rabbit mAb
Accession(Primary)	Q9BUJ2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 96 kDa; Observed MW: 118 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500

Target information (Q9BUJ2)


Gene Name	HNRNPUL1
About Gene Symbol	HNRNPUL1 (Heterogeneous Nuclear Ribonucleoprotein U Like 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of E1B AP5 (green) in Jurkat using E1B AP5 antibody, and DAPI (blue).



Western blot analysis of E1BAP5 in A549, HL-60, U251 lysates using E1BAP5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using E1BAP5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.