

## [KO]Bcr Rabbit mAb


ABZ4802

### Product Overview

|                             |                                                           |
|-----------------------------|-----------------------------------------------------------|
| <b>Name</b>                 | [KO]Bcr Rabbit mAb                                        |
| <b>Accession(Primary)</b>   | P11274                                                    |
| <b>Host Species</b>         | Rabbit                                                    |
| <b>Clonality</b>            | Monoclonal                                                |
| <b>Isotype</b>              | IgG                                                       |
| <b>Species Reactivity</b>   | Human                                                     |
| <b>Calculated MW</b>        | Calculated MW: 143 kDa; Observed MW: 160 kDa              |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:50-1:100;<br>ICC/IF-1:100-1:200 |

### Target information (P11274)

|                          |                                                                                                                                                                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | BCR                                                                                                                                                                                                                                                   |
| <b>About Gene Symbol</b> | BCR (BCR Activator Of RhoGEF And GTPase) is a GTPase-activating protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Immunohistochemistry analysis of paraffin-embedded Human Brain using Bcr antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.