

# [KO]Clathrin Heavy Chain Rabbit mAb


ABZ4803

## Product Overview


|                             |                                                |
|-----------------------------|------------------------------------------------|
| <b>Name</b>                 | [KO]Clathrin Heavy Chain Rabbit mAb            |
| <b>Accession(Primary)</b>   | Q00610                                         |
| <b>Host Species</b>         | Rabbit                                         |
| <b>Clonality</b>            | Monoclonal                                     |
| <b>Isotype</b>              | IgG                                            |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                              |
| <b>Calculated MW</b>        | Calculated MW: 192 kDa; Observed MW: 192 kDa   |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200 |

## Target information (Q00610)

|                          |                                                                                                                                                                                                                                                                                                              |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CLTC                                                                                                                                                                                                                                                                                                         |
| <b>About Gene Symbol</b> | CLTC (Clathrin Heavy Chain) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Clathrin heavy chain in HeLa, A549, HL-60, U251, U87-MG lysates using Clathrin Heavy Chain antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Clathrin heavy chain antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.