

[KD]DGKA Rabbit mAb


ABZ4804

Product Overview


Name	[KD]DGKA Rabbit mAb
Accession(Primary)	P23743
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 83 kDa; Observed MW: 83 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; IP-1:20-1:50

Target information (P23743)

Gene Name	DGKA
About Gene Symbol	DGKA (Diacylglycerol Kinase Alpha) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of DGKA in HeLa, A549 lysates using DGKA antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using DGKA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.