

[KD]PI 3 Kinase p150 Rabbit mAb


ABZ4805

Product Overview

Name	[KD]PI 3 Kinase p150 Rabbit mAb
Accession(Primary)	Q99570
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:153 kDa;Observed MW: 153 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q99570)

Gene Name	PIK3R4
About Gene Symbol	PIK3R4 (Phosphoinositide-3-Kinase Regulatory Subunit 4) is a signaling adaptor protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PI 3 Kinase p150 in Jurkat, Raji lysates using PI 3 Kinase p150 antibody.