

[KO]HARS (4H5) Mouse mAb


ABZ4811

Product Overview

Name	[KO]HARS (4H5) Mouse mAb
Accession(Primary)	P12081
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P12081)

Gene Name	HARS1
About Gene Symbol	HARS1 (Histidyl-TRNA Synthetase 1) is a biosynthetic enzyme involved in adaptive immunity. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HARS in T47D, 3T3, COS7, C6 and HeLa lysates using HARS antibody.