

[KD]Bid (3C5) Mouse mAb


ABZ4812

Product Overview


Name	[KD]Bid (3C5) Mouse mAb
Accession(Primary)	P55957
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; FC-1:50-1:100

Target information (P55957)

Gene Name	BID
About Gene Symbol	BID (BH3 Interacting Domain Death Agonist) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Bid in 293T, Jurkat and A431 lysates using Bid antibody.



Immunocytochemistry analysis of Bid (3C5) in HeLa using Bid antibody.