

[KO]ATIC Mouse mAb


ABZ4814

Product Overview

Name	[KO]ATIC Mouse mAb
Accession(Primary)	P31939
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P31939)

Gene Name	ATIC
About Gene Symbol	ATIC is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ATIC in Jurkat, C6, 3T3, COS7 and HeLa lysates using ATIC antibody.