

[KD]DDX3 Mouse mAb


ABZ4815

Product Overview

Name	[KD]DDX3 Mouse mAb
Accession(Primary)	O00571
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (O00571)

Gene Name	DDX3X
About Gene Symbol	DDX3X (DEAD-Box Helicase 3 X-Linked) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of DDX3 in COS7, Jurkat, C6, 3T3 and HeLa lysates using DDX3 antibody.