

[KO]VAMP8/EDB Rabbit mAb


ABZ4816

Product Overview

Name	[KO]VAMP8/EDB Rabbit mAb
Accession(Primary)	Q9BV40
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50


Target information (Q9BV40)

Gene Name	VAMP8
About Gene Symbol	VAMP8 (Vesicle Associated Membrane Protein 8) is a membrane receptor involved in vesicular trafficking. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of VAMP8 in HeLa, RAW264.7, NIH/3T3 lysates using VAMP8 antibody.



Western blot analysis of VAMP8 in (1) HeLa lysates; (2) NIH/3T3 lysates; (3) PC-12 lysates using VAMP8 antibody.