

[KD]Phospho-Vimentin (Ser56) Rabbit mAb


ABZ4820

Product Overview


Name	[KD]Phospho-Vimentin (Ser56) Rabbit mAb
Accession(Primary)	P08670
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:500-1:1000; FC-1:200-1:500; IP-1:20-1:50

Target information (P08670)

Gene Name	VIM
About Gene Symbol	VIM (Vimentin) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Phospho-Vimentin (Ser56) (green) in Hela using Phospho-Vimentin (Ser56) antibody, and DAPI (blue).



Western blot analysis of Phospho-Vimentin (Ser56) in U251MG lysates using Phospho-Vimentin (Ser56) antibody.