

[KD]Ephrin Receptor B3 Rabbit mAb

ABZ4823


Product Overview

Name	[KD]Ephrin Receptor B3 Rabbit mAb
Accession(Primary)	P54753
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 110 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:100


Target information (P54753)

Gene Name EPHB3

About Gene Symbol



Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of Eph receptor B3 in MOLT4 lysates using Ephrin Receptor B3 antibody.