

[KD]CDKN2A/p16INK4a Rabbit mAb


ABZ4824

Product Overview

Name	[KD]CDKN2A/p16INK4a Rabbit mAb
Accession(Primary)	P42771
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (P42771)

Gene Name	CDKN2A
About Gene Symbol	CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) is a tumor suppressor protein involved in cell cycle regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CDKN2A/p16INK4a in 293T, HeLa lysates using CDKN2A antibody.