

[KD]Carbonic Anhydrase 2 Rabbit mAb


ABZ4825

Product Overview

Name	[KD]Carbonic Anhydrase 2 Rabbit mAb
Accession(Primary)	P00918
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P00918)

Gene Name	CA2
About Gene Symbol	CA2 (Carbonic Anhydrase 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of Carbonic Anhydrase 2 in Caco 2, 293T lysates using Carbonic Anhydrase 2 antibody.



Western blot analysis of Carbonic anhydrase 2 in A431 lysates using Carbonic Anhydrase 2 antibody.