

Junctional Adhesion Molecule 1/JAM-A Rabbit mAb


ABZ4826

Product Overview


Name	Junctional Adhesion Molecule 1/JAM-A Rabbit mAb
Accession(Primary)	Q9Y624
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (Q9Y624)

Gene Name	F11R
About Gene Symbol	F11R (Junction adhesion molecule 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of JAM1 in HeLa lysates using Junctional Adhesion Molecule 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human bladder cancer using JAM1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.