

AMPK beta 1 Rabbit mAb


ABZ4827

Product Overview

Name	AMPK beta 1 Rabbit mAb
Accession(Primary)	Q9Y478
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:200-1:1000; IP-1:20-1:50

Target information (Q9Y478)

Gene Name	PRKAB1
About Gene Symbol	



Western blot analysis of AMPKβ1 in HeLa, PC-12, 3T3, COS7 lysates using AMPK beta 1 antibody.