

Axin 2 Rabbit mAb


ABZ4828

Product Overview

Name	Axin 2 Rabbit mAb
Accession(Primary)	Q9Y2T1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 94 kDa; Observed MW: 94 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100

Target information (Q9Y2T1)

Gene Name	AXIN2
About Gene Symbol	



Western blot analysis of Axin2 in (1) SW480 lysates; (2) PC-12 lysates using Axin 2 antibody.