

ATF5 Rabbit mAb


ABZ4829

Product Overview


Name	ATF5 Rabbit mAb
Accession(Primary)	Q9Y2D1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q9Y2D1)

Gene Name	ATF5
About Gene Symbol	



Western blot analysis of ATF5 in (1) Jurkat lysates; (2) 3T3 lysates; (3) C6 lysates using ATF5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast using ATF5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.