

LEF1 Rabbit mAb


ABZ4831

Product Overview


Name	LEF1 Rabbit mAb
Accession(Primary)	Q9UJU2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:100; FC-1:20; IP-1:20

Target information (Q9UJU2)

Gene Name	LEF1
About Gene Symbol	



Western blot analysis of LEF1 in Jurkat lysates using LEF1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human spleen using LEF1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.