

LEF1 Rabbit mAb


ABZ4832

Product Overview



Name	LEF1 Rabbit mAb
Accession(Primary)	Q9UJU2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 25-58 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:100; FC-1:20; IP-1:20

Target information (Q9UJU2)


Gene Name	LEF1
About Gene Symbol	




Immunocytochemistry analysis of LEF1 (green) in Jurkat using LEF1 antibody, and DAPI (blue).



Western blot analysis of LEF1 in C6 lysates using LEF1 antibody.



Western blot analysis of LEF1 in Jurkat lysates using LEF1 antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using LEF1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.