

IRAK4 Rabbit mAb


ABZ4837

Product Overview

Name	IRAK4 Rabbit mAb
Accession(Primary)	Q9NWZ3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:200; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (Q9NWZ3)

Gene Name	IRAK4
About Gene Symbol	



Western blot analysis of IRAK4 in Jurkat lysates
using IRAK4 antibody.