

ATG3 Rabbit mAb


ABZ4838

Product Overview

Name	ATG3 Rabbit mAb
Accession(Primary)	Q9NT62
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:200-1:2000

Target information (Q9NT62)

Gene Name	ATG3
About Gene Symbol	



Western blot analysis of Apg3 (Atg3) in Jurkat lysates using ATG3 antibody.