

CDT1/DUP Rabbit mAb


ABZ4844

Product Overview


Name	CDT1/DUP Rabbit mAb
Accession(Primary)	Q9H211
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 60 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (Q9H211)

Gene Name	CDT1
About Gene Symbol	CDT1 (Chromatin Licensing And DNA Replication Factor 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of CDT1 (green) in 293T using CDT1 antibody, and DAPI (blue).



Western blot analysis of CDT1 in K562, Jurkat lysates using CDT1 antibody.