

KIFC1 Rabbit mAb


ABZ4847

Product Overview


Name	KIFC1 Rabbit mAb
Accession(Primary)	Q9BW19
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:10-1:100

Target information (Q9BW19)

Gene Name	KIFC1
About Gene Symbol	



Western blot analysis of KIFC1 in (1) HeLa lysates;
(2) HepG2 lysates using KIFC1 antibody.



Immunohistochemistry analysis of paraffin-
embedded Human thyroid carcinoma using KIFC1
antibody. High-pressure and temperature Sodium
Citrate pH 6.0 was used for antigen retrieval.