

GDF15 Rabbit mAb

ABZ4848


Product Overview

Name	GDF15 Rabbit mAb
Accession(Primary)	Q99988
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 34 kDa; Observed MW: 35,13 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20

Target information (Q99988)

Gene Name GDF15

About Gene Symbol



Western blot analysis of GDF15/Mic1 in HeLa, MCF-7 lysates using GDF15 antibody.

