

## ASK1 Rabbit mAb


ABZ4850

### Product Overview


|                             |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| <b>Name</b>                 | ASK1 Rabbit mAb                                                        |
| <b>Accession(Primary)</b>   | Q99683                                                                 |
| <b>Host Species</b>         | Rabbit                                                                 |
| <b>Clonality</b>            | Monoclonal                                                             |
| <b>Isotype</b>              | IgG                                                                    |
| <b>Species Reactivity</b>   | Human, Mouse                                                           |
| <b>Calculated MW</b>        | Calculated MW: 155 kDa; Observed MW: 155 kDa                           |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:20-1:50;<br>ICC/IF-1:100-1:200; FC-1:50-1:100 |

### Target information (Q99683)

|                          |        |
|--------------------------|--------|
| <b>Gene Name</b>         | MAP3K5 |
| <b>About Gene Symbol</b> |        |



Western blot analysis of ASK1 in HeLa lysates using ASK1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using ASK1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.